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### Users are required to follow the safety and security advice.

Lubricating greases consist primarily of a base oil and a thickener. Not every base oil and not every thickener can be mixed with each other. When combining immiscible components, you get an inhomogeneous mixture that may have a low dropping point or may tend to harden. Thus, mixtures of different lubricating greases should be avoided. It is highly recommended that the original lubricating greases are used for relubricating.

The **spindle** is lubricated by us with the **OWIS® SST.F11** grease.  
Grease SST.F11 in a 5 ml applicator, part no. 90.999.0011

Relubrication should take place after about 1.000 operation hours, at least every 6 months.

The guides are maintenance free for 5 years.  
Afterwards we suggest to send the product to OWIS® for an all-encompassing maintenance.

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### Lubrication procedure

1. Move the table plate into the lower end position.
2. Turn off the control unit and make sure that it cannot be turned on accidentally.
3. Screw off the table plate (1).
4. Screw off the side panel (2).
5. Apply evenly a thin film of lubricating grease SST.F11 onto the accessible spindle surface (3). Repeat segmental by moving the stage into the upper end position until the spindle is completed overall (4). (For the complete spindle approximate 0.1 ml grease is required.)
6. Screw on the side panel, than the table plate.
7. Move the stage a few (about 5-10) times over the complete travel range to disperse the grease evenly.