


Stromlaufpläne
Pinout

SLP01

9013.0124 / 31.05.2022



44-pol HD männlich

LINPOS S und LINPOS M

	Funktion/function	Pin/Pin
Motor	Phase U	16
	Phase V	17
	Phase W	31
Hall-Sensoren/Hall sensors	5V	2
	GND	3
	Hall A+	4
	Hall A-	5
	Hall B+	6
	Hall B-	7
	Hall C+	8
	Hall C-	9
	Linearmesssystem+Endschalter/ linear measuring system + limit switches	Encoder A+
Encoder A-		11
Encoder B+		12
Encoder B-		13
Encoder I+		14
Encoder I-		15
5V		27
GND		42
MINSTOP		40
MAXSTOP		44
OWIS Steuerung/OWIS controller	OWIS ID	26
	Motorcodierung/ motor coding	28